**ATENÇÃO:**

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**ATTENTION:**

**The Technew Lapidated Crystals** is a fashion product. It's not possible to establish a precise amount of time for the crystal to remain on the teeth. Using our adhesive this period is, as a rule, 1 to 3 weeks. Using adhesive for brackets and the acid etch technique, it's possible to increase the period of permanen- ce, which will always be uncertain anyway, and which will depend on several factors, such as: some kind of hard food bitten or the collocation of the crystal closer to the incisal, etc. In case the crystal falls down the adhesive may have broken next to the tooth or next to the opaque side of the crystal. This is a random fact and does not constitute imperfection in relation to the product. The use of the same crystal for more than 3 months might cause its darkening, due to the patient's nourishing habits, such as drinking a lot of coffee, etc. In such case, the patient should be warned of the risk.

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## **Lapidated Crystals for Dental Aesthetic Technew**

### **Instructions:**

- I** - Do dental polishing with appropriate pastes;
- II** - Isolate the tooth with cotton rolls;
- III** - Acid etch the area in the tooth where the crystal will be applied;
- IV**- Apply a small amount of the adhesive of your ( the professional ) choice on the metallic base of the crystal (it is recommended to use single component adhesive or flow resin);
- V** - Put the crystal on the tooth, in the right position, let the adhesive slightly flow over the edges of the crystal in order to achieve better retention;
- VI** -Light cure it on all sides (In case a chemical curing adhesive is being used, wait until it is set.)

### **Observations:**

- I** - The Technew adhesive is used for bonding the crystal to the teeth for a short period of time (1 or 2 weeks);
- II** - In order to achieve a longer or permanent bonding between the crystal and the tooth, it is recommended that conventional bonding materials be used-acid etching the tooth, and applying the bonding system of your (the professional) choice.
- III**-Flow resin may also be used. It is recommended to let the adhesive or resin flow over the edges of the crystal in order to achieve better retention.
- IV** -The isolation of the tooth with cotton rolls should be maintained up to the end of the procedure for applying the crystal.