# Pasta Zinco-Enólica

## Indication, Use and Application of the Product

Indicated for impressions in areas without retention, especially mucosa, as for example, impression in patients without teeth.

## Specifications and technical characteristics of the product

The Zinc-enolic Paste Technew presents the following characteristics: high fluidity, small quantity of eugenol, orange flavor, low contraction during set and excellent flow. This way, the benefits of a precision impression are obtained, without tissue compression, which is ideal for patients without teeth.

The Zinc-enolic Paste Technew is an impression material obtained from the reaction of 2 pastes: the first one based on zinc oxide and the second based on Eugenol. This reaction produces the formation of zinc eugenolate, which gives the product its final characteristics, ideal for its intended use (impression material).

It is indicated for functional impressions on mouths without teeth, achieving exactness in minimum details and excellent dimensional stability. In impressions for total prosthesis, it is the election material.

Its formulation with a small amount of eugenol (less than 7%) makes it less irritant to the oral tissues. Its smooth and soft consistency allows copying the smallest details without dislocating or compressing the mucosa, allowing best adjustment of prosthesis.

## Valid date

3 years

## Use instructions

### Impression:

Dispense equal lengths of the two pastes over a glass plate or a mixture pad. That corresponds, in weight, to two parts of base paste and one part of accelerator paste.

Mix the two pastes until a uniform color of mixture is obtained. It takes around 1 minute.

Put the mixture inside a prepared individual tray, which must be dry.

With a gentle pressure, put the tray inside the mouth. Register the muscular insertions for about 30 to 45 seconds after inserting the tray inside the mouth.

Obs.: The oral tissue must be clean and salivation controlled. Lubricate the lips and adjacent areas in order not to allow the mixture to stick to them.

Keep a gentle pressure until the mixture sets (around minutes after starting mixing), and then remove the impression.

## Work time: 1 minute

Set time: 3 minutes

## Dropping the stone:

Wash the impression obtained with cool water to remove oral fluids. Take away the excess water by agitating or using an air spray.

Mix the stone and pour it over the impression with a short vibration. It is not necessary to use a separating medium.

After the final set of the stone, immerse the set in water enough warmed to allow the impression to be free from the retentions. Do not let water boil.

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